

# Evaluation of the R-Biopharm RIDA® GENE PVL real-time PCR assay for the laboratory diagnosis of skin and soft tissue infections caused by PVL

Andreas Hiergeist<sup>1</sup>, Andreas Simons<sup>2</sup>,  
<sup>1</sup> microBIOMix GmbH, Regensburg, Germany  
<sup>2</sup> R-Biopharm AG, Darmstadt, Germany

## Objectives

*Staphylococcus aureus* (SA) is a major cause of nosocomial infections in hospitals and healthcare settings. Since the mid-1990s the number of infections in the population increased with no previous history of medical facility contact. This increase in infections in the population is caused by *Staphylococcus aureus* strains that carry the virulence factor Pantone-Valentine leukocidin (PVL). PVL is a cytotoxin encoded by the lukF-PV and lukS-PV genes and contributes to tissue necrosis by lysis of macrophages and neutrophils hence inhibiting the process of wound closure. PVL can be produced by methicillin-sensitive *Staphylococcus aureus* strains (MSSA) as well as methicillin-resistant *Staphylococcus aureus* strains (MRSA). MRSA strains that carry the virulence factor PVL are called community-acquired (CA)-MRSA.<sup>1,2</sup> Risk groups for transmission of CA-MRSA or PVL-MSSA are for example families, persons performing close contact sports, persons from educational settings, prisoners and military personnel.<sup>3,4</sup> The clinical manifestation of PVL-positive *Staphylococcus aureus* (PVL-SA) strains are skin and soft tissue infections, particularly recurrent invasive abscesses and in the US, CA-MRSA is already the most common pathogen causing this clinical picture in outpatients.<sup>2</sup> In the diagnosis of skin and soft tissue infections, but also in systemic diseases caused by PVL-SA, it is useful to determine the virulence factor PVL to initiate appropriate therapeutic and hygienic measures. This study evaluated the performance of the RIDA® GENE PVL assay for the detection of skin and soft tissue infections caused by PVL compared to an established in-house real-time PCR assay.

## Methods

To evaluate the sensitivity and specificity of the RIDA® GENE PVL real-time PCR assay, 20 known *Staphylococcus aureus*-positive nasal swabs, quantitatively spiked with DNA of PVL-positive strains and 30 known *Staphylococcus aureus*-negative nasal swabs were tested. In a parallel study, 32 known PVL-SA samples (17 MRSA-SA and 15 MSSA-SA samples) and

62 PVL-negative MRSA/MSSA samples were grown on blood agar plates (Figure 1). For both studies, an established in-house real-time PCR assay was defined as gold standard. DNA was analysed on the LightCycler® 480II (Roche). The RIDA® GENE PVL real-time PCR assay is a qualitative assay using hydrolysis probes for the qualitative detection of PVL.

Figure 1 Study design



Table 1 Analyte target gene and detection channel overview

Analyte	Target Gene	Detection Channel Light® Cycler 480II
PVL	lukF-PV	465/510
ICD	Synthetic target DNA sequence	533/580

## Results

Using the in-house real-time PCR as gold standard, all 20 known *Staphylococcus aureus*-positive samples, as well as 30 known *Staphylococcus aureus*-negative samples from nasal swabs were correctly identified by the RIDA® GENE PVL real-time PCR assay indicating a sensitivity and specificity of 100 % for the RIDA® GENE

PVL real-time PCR assay (Figure 2). Of the 94 samples from culture, 4 samples showed an inhibition in both PCR methods. Compared to the in-house real-time PCR, the RIDA® GENE PVL real-time PCR assay was in 100 % agreement (Figure 3).

Picture 1 RIDA® GENE PVL



Figure 2 In-house RIDA® GENE PVL vs. real-time PCR from nose swabs

		In-house real-time PCR		Total
		+	-	
RIDA® GENE PVL	+	20	0	20
	-	0	30	30
Total		20	30	50

Sensitivity:	100.0 %
Specificity:	100.0 %
PPV:	100.0 %
NPV:	100.0 %

Figure 3 In-house RIDA® GENE PVL vs. real-time PCR from culture

		In-house real-time PCR		Total
		+	-	
RIDA® GENE PVL	+	32	0	32
	-	0	62	62
Total		32	62	94

Sensitivity:	100.0 %
Specificity:	100.0 %
PPV:	100.0 %
NPV:	100.0 %

## References:

- David MZ and Daum RS. Community-Associated Methicillin-Resistant *Staphylococcus aureus*: Epidemiology and Clinical Consequences of an Emerging Epidemic. CLINICAL MICROBIOLOGY REVIEWS 2010, 23(3): 616–687.
- Linde HJ and lehn N. Community-associated MRSA: Klinik, Therapie, Hygiene. Krankenh.hyg. update 2008, 3(1):29-44.
- Zetola N et al. Community-acquired methicillin resistant *Staphylococcus aureus*: an emerging threat. Lancet Infect Dis 2005, 5: 275-286.

## Conclusion

The RIDA® GENE PVL assay is a highly sensitive and rapid qualitative real-time PCR assay for the detection of the PVL gene from human nasal swabs and culture.

The RIDA® GENE PVL real-time PCR shows good correlation with an established real-time PCR for the detection of *Staphylococcus aureus* of skin and soft tissue infections caused by PVL.

Results are available in less than 2 hours.