



Need for IFX and ADM rapid testing? Therapeutic drug monitoring (TDM)

The lateral flow assays RIDA®QUICK IFX Monitoring and RIDA®QUICK ADM Monitoring are rapid assays, which allow the quantitative determination of infliximab and adalimumab trough levels (TLs) within 20 minutes. The assays facilitate the use of therapeutic drug monitoring in hospitals and infusion centers and the decision making in regards to therapy adaptation.

RIDA®QUICK IFX Monitoring, Art. No. GN3041



- Fast and reliable drug monitoring of infliximab (IFX) and its biosimilars as well as adalimumab (ADM)
- Quantitative lateral flow tests
- 20 minutes time to result
- High correlation to the corresponding ELISA assays RIDASCREEN® IFX Monitoring and RIDASCREEN® ADM Monitoring due to identical monoclonal antibodies (MA-IFX6B7 and MA-ADM40D8)
- 25 test cassettes per package
- 2 controls available as accessory

RIDA®QUICK IFX Monitoring and the RIDA®QUICK ADM Monitoring are read on the RIDA®QUICK SCAN II for precise and accurate measurement.

Scan to view video clip about
RIDA®QUICK IFX Monitoring
([https://clinical.r-biopharm.com/
products/ridaquick-ifx-monitoring/](https://clinical.r-biopharm.com/products/ridaquick-ifx-monitoring/))



Products for therapeutic drug monitoring (TDM)

Lateral flow assays

Product	Description	Tests	Matrix	Art. No.
Lateral flow test and device				
RIDA®QUICK IFX Monitoring	Immunochromatographic lateral flow assay for the quantitative determination of infliximab (IFX) and its biosimilars	25	Serum/plasma	GN3041
New RIDA®QUICK ADM Monitoring	Immunochromatographic lateral flow assay for the quantitative determination of adalimumab (ADM)	25	Serum/plasma	GN3043
RIDA®QUICK IFX Monitoring Control Set	Positive controls • accessory for GN3041			GP3041
New RIDA®QUICK ADM Monitoring Control Set	Positive controls • accessory for GN3043			GP3043
RIDA®QUICK SCAN II - IVD SET	Lateral flow reader (CE-IVD) and 2D barcode scanner • for read out of GN3041 and GN3043			ZRQS2-KD-SET

ELISAs based on assays of KU Leuven, Belgium

Enzyme immunoassay				
RIDASCREEN® IFX Monitoring	Enzyme immunoassay for the quantitative determination of infliximab (IFX) and its biosimilars	96	Serum/plasma	G09041
RIDASCREEN® Anti-IFX Antibodies	Enzyme immunoassay for the quantitative determination of antibodies to infliximab (IFX) and its biosimilars	96	Serum/plasma	G09042
RIDASCREEN® ADM Monitoring	Enzyme immunoassay for the quantitative determination of adalimumab (ADM)	96	Serum/plasma	G09043
RIDASCREEN® Anti-ADM Antibodies	Enzyme immunoassay for the quantitative determination of antibodies to adalimumab (ADM)	96	Serum/plasma	G09044
RIDASCREEN® VDZ Monitoring	Enzyme immunoassay for the quantitative determination of vedolizumab (VDZ)	96	Serum/plasma	G09045
New RIDASCREEN® GLM Monitoring	Enzyme immunoassay for the quantitative determination of golimumab (GLM)	96	Serum/plasma	G09047

For further details and information visit our website www.r-biopharm.com, contact your local distributor or Clinical Sales International.

R-Biopharm contacts:

Clinical Sales International:
Phone: +49 (0) 61 51 - 81 02-0
Fax: +49 (0) 61 51 - 81 02-688
E-mail: clinical.sales@r-biopharm.de

Order Department:
Phone: +49 (0) 61 51 - 81 02-0
Fax: +49 (0) 61 51 - 81 02-20
E-mail: orders@r-biopharm.de