



# RIDA® QUICK SCAN II

Easy and quick quantitative lateral flow measurements

Our new LTF reader brings your measurements to a more precise and accurate level

- High resolution image evaluation
- Methods and test assay specifications via QR code technology
- LED self-calibration guarantees signal stability and independent measurements from outside factors like dust
- Maintenance free
- Operating time (battery): 6 h
- Up to 1400 test results can be stored
- Results as PDF, TXT and XML files
- LAN connectivity
- 4.3" color touchscreen
- Developed for R-Biopharm's Clinical Diagnostics and Food & Feed Analysis

RIDA® QUICK SCAN II loading process



## Product information

Product	Description	Packaging size	Art. No.
RIDA®QUICK SCAN II – IVD Set	Lateral flow reader (IVD, CE) and 2D barcode scanner	1	ZRQS2-KD-SET
RIDA®QUICK SCAN II	Lateral flow reader (IVD, CE)	1	ZRQS2-KD
Honeywell Xenon 1900	2D barcode scanner for ZRQS2-KD/FF	1	ZRQS2-BS

For further details and information visit our website [www.r-biopharm.com](http://www.r-biopharm.com), contact your local distributor or Clinical Sales International.

### R-Biopharm contacts:

**Clinical Sales International:**  
Phone: +49 (0) 61 51 - 81 02-0  
Fax: +49 (0) 61 51 - 81 02-688  
E-mail: [clinical.sales@r-biopharm.de](mailto:clinical.sales@r-biopharm.de)

**Order Department:**  
Phone: +49 (0) 61 51 - 81 02-0  
Fax: +49 (0) 61 51 - 81 02-20  
E-mail: [orders@r-biopharm.de](mailto:orders@r-biopharm.de)